



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Lab Coats

**Product Code :**  
LBNY-00017-11000110



**Description :**

Provide excellent protection against submicron particles.  
Lightweight and tea-resistant Tyvek fabric  
Knee length with five-snap front and full sleeves  
Disposable, yet strong enough to withstand several wearings.  
Provide limited splash protection  
Chemically and biologically inert  
White  
Applications:  
Lead and asbestos abatement/remediation, general maintenance/operations, spray painting, general clean-up.

**Technical Specification :**

Lab Coats



**Laboratory Instrument India**