



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Lab Coats

**Product Code :**  
LBNY-00017-11000110



**Description :**

Provide excellent protection against submicron particles.

Lightweight and tear-resistant Tyvek fabric

Knee length with five-snap front and full sleeves

Disposable, yet strong enough to withstand several wearings.

Provide limited splash protection

Chemically and biologically inert

White

Applications:

Lead and asbestos abatement/remediation, general maintenance/operations, spray painting, general clean-up.

**Technical Specification :**

Lab Coats



**Laboratory Instrument India**