



Email : sales@laboratoryinstrumentindia.com

Product Name :
Lab Coats

Product Code :
LBNY-00017-11000110



Description :

Provide excellent protection against submicron particles.
Lightweight and tea-resistant Tyvek fabric
Knee length with five-snap front and full sleeves
Disposable, yet strong enough to withstand several wearings.
Provide limited splash protection
Chemically and biologically inert
White
Applications:
Lead and asbestos abatement/remediation, general maintenance/operations, spray painting, general clean-up.

Technical Specification :

Lab Coats



Laboratory Instrument India