



Email : sales@laboratoryinstrumentindia.com

Product Name :
Eyes Only Cleanroom Hoods

Product Code :
LBNY-00017-11000111



Description :

Clean-processed garments offer low level of particle shedding
Cleanroom hood features bound seams for strength and durability. Available in clean or clean and sterile versions.
Eyes-only opening
Snaps for adjustable fit
Bound seams are covered with garment fabric for reinforcement
Bound hood opening for lower particle shedding
Available in clean or clean and sterile versions
Available in medium or large sizes
Greener Choice Claims
Environmental benefits include:
Recyclable.

Technical Specification :

Eyes Only Cleanroom Hoods



Laboratory Instrument India

