



Email : sales@laboratoryinstrumentindia.com

Product Name :
Liquid and Particle Protection Boot Covers

Product Code :
LBNY-00017-11000112



Description :

Elastic opening at top and elastic ankle
Serged seams
Highly durable PVC sole
Made to keep out debris, dry particulates, and liquid splashes
Boot covers feature a highly durable PVC sole and are made with breathable, microporous film to keep out debris, non-hazardous liquid chemicals, and dry particulates.
Boot Covers available in 3 sizes, choose from 17 in. or 18 in. height
Adheres to criteria for antistatic materials.

Technical Specification :

Liquid and Particle Protection Boot Covers



Laboratory Instrument India