



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Frock with Mandarin Collar

**Product Code :**  
LBNY-00017-11000114



**Description :**

Ideal for lead and asbestos abatement/remediation, general maintenance/operations, spray painting, general clean-up.

Flash-spun high-density polyethylene provides an ideal balance of protection, durability and comfort

Body-length jacket with five front snap closures for easy donning and doffing

Mandarin collar fits tightly around hooded PAPR

Elastic opening for tighter fit at wrist

Enhanced frontal protection

Inherent barrier to small size particles built into the fabric itself — no films or laminates to abrade or wear away

White

Serged seams.

**Technical Specification :**

Frock with Mandarin Collar



**Laboratory Instrument India**

