



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Vinyl Aprons, Blue, Sewn Edges

**Product Code :**  
LBNY-00017-11000115



**Description :**

Suitable when lighter aprons will do  
Grommets and bib are electronically heat-sealed for added strength  
Nylon ties secure apron at back  
Available with nylon-serged finished edges.

**Technical Specification :**

Vinyl Aprons, Blue, Sewn Edges



**Laboratory Instrument India**