



Email : sales@laboratoryinstrumentindia.com

Product Name :
Nonsterile Sleeves

Product Code :
LBNY-00017-11000117



Description :

Comfortable, lightweight, and durable Tyvek sleeves are available in bulk or clean packaging.
Comfortable, lightweight and durable.
Covered elastic wrist, bicep
Bound seams
White
Length: 18 in. (45.7cm)
Greener Choice Claims
Environmental benefits include:
Recyclable.

Technical Specification :

Nonsterile Sleeves



Laboratory Instrument India