



Email : sales@laboratoryinstrumentindia.com

Product Name :
Serged Seam Cleanroom Frocks

Product Code :
LBNY-00017-11000118



Description :

Comfortable, lightweight, and durable
Bulk packaged in double transparent poly liners
Snap closure (6 snaps)
Elastic wrists
Excellent barrier to dry particulates, nonhazardous liquids, and microorganisms.
No pockets to trap contamination
Color: White
Greener Choice Claims
Environmental benefits include:
Recyclable.

Technical Specification :

Serged Seam Cleanroom Frocks



Laboratory Instrument India

