



Email : sales@laboratoryinstrumentindia.com

Product Name :
Coveralls, Attached Hood

Product Code :
LBNY-00017-11000120



Description :

Use in controlled environments with high standards for particle and microbiological contamination control
Pharmaceutical, biotech, medical device, electronic settings
Protects against micro-organisms and particulates.
Spinning process creates barrier to particles, micro-organisms, and nonhazardous liquids
Light and durable coveralls provide a barrier to particles, microorganisms, and non-hazardous light liquid splash.
Elastic wrists and ankles
Zipper closure
Lowest level particle shedding
Standard hood
Recyclable
Serged seams
Greener Choice Claims
Environmental benefits include:
Recyclable.

Technical Specification :

Coveralls, Attached Hood

Laboratory Instrument India

