




Email : sales@laboratoryinstrumentindia.com

Product Name :
Breathable Particle Protection Bouffant Cap

Product Code :
LBNY-00017-11000122



Description :

Made with breathable Microforce barrier fabric
Features serged seams and an elastic opening
Lightweight and breathable 24 in. bouffant caps are best suited for use in cleanroom applications.
Recyclable via The Right Cycle Program by Kimberly-Clark Professional
Elastic is not made with natural rubber latex (reducing risk for individuals with a Type 1 latex allergy)
Passes for antistatic (anti stat) materials
Greener Choice Claims
Environmental benefits include:
Recyclable.

Technical Specification :

Breathable Particle Protection Bouffant Cap



Laboratory Instrument India