



Email : sales@laboratoryinstrumentindia.com

Product Name :
Frock with Mandarin Collar

Product Code :
LBNY-00017-11000124



Description :

Flash-spun high-density polyethylene provides an ideal balance of protection, durability and comfort
Body-length jacket with five front snap closures for easy donning and doffing
Enhanced frontal protection
Inherent barrier to small size particles built into the fabric itself — no films or laminates to abrade or wear away
Ideal for lead and asbestos abatement/remediation, general maintenance/operations, spray painting, general clean-up.
Elastic opening for tighter fit at wrist
Mandarin collar fits tightly around hooded PAPR
Serged seams
White.

Technical Specification :

Frock with Mandarin Collar



Laboratory Instrument India

