



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Frock with Mandarin Collar

**Product Code :**  
LBNY-00017-11000124



**Description :**

Flash-spun high-density polyethylene provides an ideal balance of protection, durability and comfort  
Body-length jacket with five front snap closures for easy donning and doffing  
Enhanced frontal protection  
Inherent barrier to small size particles built into the fabric itself — no films or laminates to abrade or wear away  
Ideal for lead and asbestos abatement/remediation, general maintenance/operations, spray painting, general clean-up.  
Elastic opening for tighter fit at wrist  
Mandarin collar fits tightly around hooded PAPR  
Serged seams  
White.

**Technical Specification :**

Frock with Mandarin Collar



**Laboratory Instrument India**

