



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
PH Meter

**Product Code :**  
LBNY-0001-980004



**Description :**

PH Meter

**Technical Specification :**

**Fast and Accurate:** No need to work with color charts or messy test strips, this meter get the result fast and accurate. It with 0.01 resolution, 0.05 accuracy, full pH 0-14 measurements! And this ph meter is factory calibrated. And when the product with problems, users can make calibration to get problem solved.

**Wide Use:** This PH meter is great for a wide range.

**Easy calibration:** No need any adjustment screwdrivers, just press a button to finish the calibration.

**ATC Function:** ATC accounts for changes in temperature automatically, it can make the meter more accurate.



**Laboratory Instrument India**