



Email : sales@laboratoryinstrumentindia.com

Product Name :
Long Sleeve Aprons

Product Code :
LBNY-00017-11000132



Description :

Provide protection from a wide variety of chemical environments

Frontal protection

Protection in chemical mixing, hazardous remediation, emergency medical response, paint spraying and nuclear environments.

30 minutes of protection in 124 chemical situations

Film laminated to Tyvek resists damage

Open back, elastic wrist, and reinforced seams

Ties at waist, elastic waist opening.

Technical Specification :

Long Sleeve Aprons



Laboratory Instrument India