



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sterile Coveralls with Hood and Thumb Loops

Product Code :
LBNY-00017-11000133



www.laboratoryinstrumentindia.com

Description :

Disposable sterile coverall with hood features a thumb loop for secure hold.
Elasticated back, cuffs and ankles
Anti-Static properties
Low linting, anti-static Fabric made
Sterilized by Gamma irradiation
Zip front with flap cover
Thumb loops on wrist
Low-particulate
Sterility Assurance Level 10-6
Manufactured in a facility holding
Available in sizes ranging from S through 4XL.

Technical Specification :

Sterile Coveralls with Hood and Thumb Loops



Laboratory Instrument India

