



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Clean Processed Cleanroom Coveralls, Individually  
Packaged

**Product Code :**  
LBNY-00017-11000135



**Description :**

Coveralls feature a polypropylene breathable SMS fabric with a cloth-like feel.  
The Reflex Coverall Design is 7.5 times more rip-resistant than ANSI minimum standards and offers excellent fit and freedom of movement for the wearer.  
This personal protective apparel is triple-bagged and vacuum packaged, suitable for ISO Class 5 or higher cleanrooms.  
Clean Processed Cleanroom Coveralls/Case  
Triple bagged, vacuum packed, suitable for cleanrooms  
Reflex Coverall Design offers excellent fit and freedom of movement and 7.5 times more rip-resistant than ANSI minimums  
Clean Processed Cleanroom Coveralls have an average bacterial filtration efficiency (BFE) of 96% and are made of 100% polypropylene breathable SMS fabric, which gives the wearer cloth-like comfort.  
It is made without latex, silicone and the BHT preservative and is designed for single-use (disposable PPE). The easy-to-don, individually packaged coveralls feature a mandarin collar, thumb loops and tunneled elastics.  
Clean-processes and individually packaged under cleanroom conditions  
This protective apparel isn't made with natural rubber latex, BHT preservative or silicone  
Greener Choice Claims  
Environmental benefits include:  
Recyclable.

**Technical Specification :**

Clean Processed Cleanroom Coveralls, Individually Packaged



**Laboratory Instrument India**