



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Disposable Long Sleeve Cover

**Product Code :**  
LBNY-00017-11000137



**Description :**

Extra long length sleeve covers are constructed from antistatic, low-linting Low linting, anti-static Clean Tough fabric  
Elasticated arms and cuffs  
Clean Tough material and are ISO Class 4 compatible.  
Protective sleeves available in sterile or non-sterile  
Universal size  
Overall length: 500 mm

**Technical Specification :**

Disposable Long Sleeve Cover



**Laboratory Instrument India**