



Email : sales@laboratoryinstrumentindia.com

Product Name :
Disposable Long Sleeve Cover

Product Code :
LBNY-00017-11000137



Description :

Extra long length sleeve covers are constructed from antistatic, low-linting Low linting, anti-static Clean Tough fabric
Elasticated arms and cuffs
Clean Tough material and are ISO Class 4 compatible.
Protective sleeves available in sterile or non-sterile
Universal size
Overall length: 500 mm

Technical Specification :

Disposable Long Sleeve Cover



Laboratory Instrument India