



Email : sales@laboratoryinstrumentindia.com

Product Name :
127 Coveralls

Product Code :
LBNY-00017-11000138



Description :

Flexible and durable light liquid splash protection
Elastic opening for tighter fit at ankle
Elastic opening for tighter fit at wrist
Attached hood with elastic around face opening
Yellow for high visibility
Polyethylene-coated fabric provides at least 30 minutes of protection against 42 chemical challenges
10 mil thickness.

Technical Specification :

127 Coveralls



Laboratory Instrument India