



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cleanroom Non-Sterile Sleeve Protectors

Product Code :
LBNY-00017-11000145



Description :

Made with advanced barrier fabric to provide enhanced protection against liquid exposure
Double bagged in bulk packaging for convenience and waste reduction
Recyclable through the RightCycle Program by Professional
Cleanroom sleeve protectors are made with advanced barrier fabric to provide enhanced protection against liquid exposure in cleanroom environments.
Reinforced seams and thumb loops for full wrist coverage
Folded inside out for aseptic donning, including a blue indicator line on the interior
200 Cleanroom Sleeve Protectors/Case, 4 Bags of 50 Sleeves

Technical Specification :

Cleanroom Non-Sterile Sleeve Protectors



Laboratory Instrument India