



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Bib Apron

**Product Code :**  
LBNY-00017-11000146



**Description :**

The inherent barrier stops submicron particles even after abrasion.  
Provides enhanced frontal protection  
Nonwoven Tyvek fabric helps resist penetration by particles and other contaminants.  
Bound seams around neck  
One size fits most  
Ties at waist

**Technical Specification :**

Bib Apron



**Laboratory Instrument India**