



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Bound Seam Coveralls, Clean-Processed

**Product Code :**  
LBNY-00017-11000148



**Description :**

Clean-processed coveralls provide an inherent barrier to particles, microorganisms, and non-hazardous light liquid splash. Bound seams reduce the potential for particle penetration.  
Elastic openings for tighter fit at wrist and ankle  
Clean-processed garments offer lowest level of particle shedding.  
Bound seams are covered with garment fabric to reinforce the seam  
Sleeve design for greater range of motion and comfort  
Bound neck for lower particle shedding  
Front zipper closure for easy donning and doffing.

**Technical Specification :**

Bound Seam Coveralls, Clean-Processed



**Laboratory Instrument India**