



Email : sales@laboratoryinstrumentindia.com

Product Name :

Breathable Particle Protection Coveralls: Elastic Back, Wrist and Ankles with Attached Hood and Boots

Product Code :

LBNY-00017-11000150



Description :

Air and moisture vapor can pass through to aid wearer comfort

Triple-stitched serged seams

Patented Reflex design is 7.5 times less likely to rip out than minimums — provides 12% more chest room and 6% longer body length

Provide general protection against many particles

Breathable patented barrier SMS fabric passes NFPA 99 for antistatic materials

For use in laboratories, remediation, decontamination modesty garment, pesticides, and heavy dust environments. Elastic Back, Wrist and Ankles with Attached Hood and Boots.

Offers high abrasion resistance and is low linting

Withstand multiple wearings

White.

Technical Specification :

Breathable Particle Protection Coveralls: Elastic Back, Wrist and Ankles with Attached Hood and Boots



Laboratory Instrument India

