



Email : sales@laboratoryinstrumentindia.com

Product Name :
Shoe Covers, Bulk

Product Code :
LBNY-00017-11000153



Description :

Flash spun polyethylene creates barrier to particles, micro-organisms and non-hazardous liquids
Elastic opening
Serged seams
Suitable for use in controlled environments with high standards for particle and microbiological contamination control (pharmaceutical, biotech, medical device and electronic settings).
Gripper PVC soles
Should not be worn around open flame, heat source, sparks or ignition sources.

Technical Specification :

Shoe Covers, Bulk



Laboratory Instrument India