



Email : sales@laboratoryinstrumentindia.com

Product Name :
Full Face Opening Hoods

Product Code :
LBNY-00017-11000155



Description :

Full face opening
Ties with loops for adjustable fit
Bound seams are covered with garment fabric for reinforcement
Bound hood opening for lower particle shedding
Open face hoods provide an inherent barrier to particles, microorganisms, and non-hazardous light liquid splash.
Available in non-sterile, clean-processed, and sterile versions.
One size fits most
Available in bulk packaging with double transparent poly liners
Available in non-sterile, clean-processed, and sterile versions
Clean-processed garments Class Cleanroom 4 standards.

Technical Specification :

Full Face Opening Hoods



Laboratory Instrument India

