



Email : sales@laboratoryinstrumentindia.com

Product Name :
Thermal Conductor

Product Code :
LBNY-0001-980007



Description :

Thermal Conductor

Technical Specification :

The apparatus is designed to determine thermal conductivity of different liquids. The apparatus consists of a heater. The heater heats a thin layer of liquid. A cooling plate removes heat through liquid layer, ensuring unidirectional heat flow. Temperature is measured across the liquid layer and complete assembly is properly insulated. A proper arrangement for changing the liquids is provided. The whole assembly is kept in a chamber.



Laboratory Instrument India