



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
QC Sleeve

**Product Code :**  
LBNY-00017-11000156



**Description :**

A lightweight, and durable fabric, utilizes the strength of DuPont Tyvek fabric and a polyethylene coating. fabric provides at least 30 minutes of protection against 42 chemical challenges.

Sleeve. 18 in. Long. Elastic at Both Ends. Bound Seams. Yellow.

Used for light splash protection in a variety of industrial environments, including petroleum refining, pulp and paper manufacturing, food processing, chemical processing, and pharmaceutical manufacturing.

**Technical Specification :**

QC Sleeve



**Laboratory Instrument India**