



Email : sales@laboratoryinstrumentindia.com

Product Name :
Coveralls with Hood and Boot

Product Code :
LBNY-00017-11000159



Description :

Made with a microporous film that is laminated to a spund bonded polypropylene substrate
Standard hood with elastic around the face
Made from microporous film, 50 garments help protect workers who face on-the-job exposure to a range of non-hazardous aerosols, light liquid splash and dry particles.
Elastic wrist
Serged seams
Attached skid-resistant boots
Front zipper closure
White.

Technical Specification :

Coveralls with Hood and Boot



Laboratory Instrument India