



Email : sales@laboratoryinstrumentindia.com

Product Name :
30 Shoe Covers

Product Code :
LBNY-00017-11000160



Description :

Material polyethylene coated spun bond polypropylene
Generously sized in large or x-large
Resistant to liquids
Lint-free
Elastic ankles for secure fit
Will not delaminate
Choose from blue or white slip-resistant sole shoe covers suitable for Cleanroom ISO 6 (Class 1000) controlled environments.
Slip resistant sole shoe covers
Available in blue or white
Seamless
Protects against non-hazardous dry particulates, liquids, and powders (Class 1,000) and above.

Technical Specification :

30 Shoe Covers

Laboratory Instrument India

