



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Breathable Apron

**Product Code :**  
LBNY-00017-11000162



**Description :**

Made from breathable, patented barrier SMS fabric.  
Secure fit  
Ties in back  
Measures 40L x 28 in. W (102 x 71cm)  
Passes NFPA 99 criteria for antistatic materials  
Available in denim blue or white.

**Technical Specification :**

Breathable Apron



**Laboratory Instrument India**