



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Pullover Hood

**Product Code :**  
LBNY-00017-11000163



**Description :**

Features tunnelized elasticized face opening and provides complete coverage of head and shoulders  
Lightweight, chemically and biologically inert and tear resistant  
Excellent for use with coveralls  
Spunbound olefin hood provides protection against microscopic airborne particles.  
One-size-fits-all  
White.

**Technical Specification :**

Pullover Hood



**Laboratory Instrument India**