



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Light Duty Shoe Covers

**Product Code :**  
LBNY-00017-11000165



**Description :**

Made of lightweight spunbond polypropylene material  
Fully elasticized and keep workers cool and comfortable  
For a wide range of work environments such as food processing, manufacturing facilities, hospitals, and construction sites.  
Universal size  
Blue.

**Technical Specification :**

Light Duty Shoe Covers



**Laboratory Instrument India**