



Email : sales@laboratoryinstrumentindia.com

Product Name :
Coveralls, Clean-Processed and Sterile

Product Code :
LBNY-00017-11000167



Description :

Coverall has elastic openings for tighter fit at wrist and ankle
Bound neck for lower particle shedding
Front zipper closure for easy donning and doffing
Sieve design for greater range of motion and comfort
Clean-processed and sterile coveralls worn in controlled environments protect against particulates and other contaminants.
Clean-processed garments offer lowest level of particle shedding.
Bound seams are covered with garment fabric to reinforce the seam and to reduce the potential for particle penetration
Garments are gamma sterilized to SAL of 10⁻⁶
Full traceability on all sterilized apparel with sterility available.

Technical Specification :

Coverall, Clean-Processed and Sterile



Laboratory Instrument India

