



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Sterile Cleanroom Sleeve Protector

**Product Code :**  
LBNY-00017-11000169



**Description :**

Made of 100% polypropylene breathable SMS fabric  
Innovative cuff fold, which presents the inside of the sleeves as the package is opened, reducing the risk of contaminating the outside of the sleeve  
Integrated thumb loop keeps the sleeve from rolling back and exposing the arm or wrist  
Sterile sleeve protector with integrated thumb loop are triple bagged, vacuum packed, and suitable for Class 5 or higher cleanrooms.  
sterile cleanroom sleeve protectors have an average bacterial filtration efficiency (BFE) of 96%  
Designed for single-use  
100 Pairs (200 each) per case  
Triple-bagged and vacuum packaged  
Not made with natural rubber latex, silicone, or the BHT preservative  
Each individually bagged pair of sleeves has a sterilization indicator as a visual that it's undergone gamma sterilization.

**Technical Specification :**

Sterile Cleanroom Sleeve Protector

---

# Laboratory Instrument India

