



Email : sales@laboratoryinstrumentindia.com

Product Name :
2-1-2 Coveralls, Sterile

Product Code :
LBNY-00017-11000174



Description :

Made from flash spun polyolefin protective material to create a barrier to particulates, nonhazardous liquids and micro-organisms

Ideal protection against micro-organisms and particulates in controlled environments.

Coveralls feature Micro-Clean with 2-1-2 polymeric resin coating that allows for thorough protection and reduction in static.

Antistatic

Bound neck and seams with zipper closure

Dolman sleeves

Elastic wrists and ankles.

Available clean processed (individually packaged) or sterilized via gamma irradiation

Coated in 2-1-2 polymeric resin.

Technical Specification :

2-1-2 Coveralls, Sterile



Laboratory Instrument India

