



Email : sales@laboratoryinstrumentindia.com

Product Name :
Boot Covers Bound

Product Code :
LBNY-00017-11000175



Description :

Made from flash spun polyolefin protective material using proprietary flash spinning process
Comfortable, lightweight and durable
Provides inherent barrier to dry particles, microorganisms, and non-hazardous liquids. Available in bulk or clean-processed packaging.
Covered elastic opening
Available with choice of seams, openings and soles
Ankle ties for secure fit
Height: 18 in. (45.7cm)
White.

Technical Specification :

Boot Covers Bound



Laboratory Instrument India