



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Boot Covers Bound

**Product Code :**  
LBNY-00017-11000175



**Description :**

Made from flash spun polyolefin protective material using proprietary flash spinning process  
Comfortable, lightweight and durable  
Provides inherent barrier to dry particles, microorganisms, and non-hazardous liquids. Available in bulk or clean-processed packaging.  
Covered elastic opening  
Available with choice of seams, openings and soles  
Ankle ties for secure fit  
Height: 18 in. (45.7cm)  
White.

**Technical Specification :**

Boot Covers Bound



**Laboratory Instrument India**