



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sterile Cleanroom Boots, Ties

Product Code :
LBNY-00017-11000178



Description :

Made without latex, silicone or the BHT preservative and are designed for single-use
Made of 100% polypropylene breathable SMS fabric
Sterile cleanroom boots with ties, vinyl
Average bacterial filtration efficiency (BFE) of 96%
Sterile cleanroom boots with anti-slip soles feature innovative CLEAN-DON technology for aseptic donning.
Suitable for EU GMP ISO 5 Grade A sterile cleanrooms.
Triple-bagged and vacuum packaged
Clean-Don ties and 2-inch edge vinyl, for an enhance moisture barrier
100 Pairs/200 Each
Cleanroom Apparel is manufactured, clean-processes and individually packaged under cleanroom conditions.
Each individually bagged pair of boot covers has a sterilization indicator as a visual that it's undergone gamma sterilization.

Technical Specification :

Sterile Cleanroom Boots, Ties



Laboratory Instrument India

