



Email : sales@laboratoryinstrumentindia.com

Product Name :
Protective Coveralls: With Elastic Wrists and Ankles

Product Code :
LBNY-00017-11000179



Description :

Disposable coveralls can be easily discarded via landfill burial or proper incineration³

Long torso with polyester front zipper measures 26 in. long (66cm)

White.

Lightweight coveralls (1.8 oz. per sq. yd.) provide excellent barrier protection against nontoxic particle, liquid and spray hazards. True aerosol- and liquid-proof microporous barrier demonstrates >99% holdout of particles down to 0.35 μ m.

Affordable high-barrier protection

Proprietary bilaminate composite fabric provides superior liquid penetration¹

Compared favorably to leading fabric in physical property tests covering resistance to ripping and tearing²

Breathable, cloth-like material treated to provide antistatic properties for use in controlled environments

Low linting for use in critical applications.

Technical Specification :

Protective Coveralls: With Elastic Wrists and Ankles



Laboratory Instrument India

