



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sterile Integrated Hood and Mask

Product Code :
LBNY-00017-11000180



Description :

Polypropylene breathable SMS fabric
Clean-don ties for easy gowning and are folded inside out for aseptic donning
Protective hood-and-mask-combos are made without latex, silicone or the BHT preservative
Average bacterial filtration efficiency (BFE) of 96%
Cleanroom sterile hood with a mask integrated into the hood.
Triple-bagged and vacuum packaged
Designed for single-use
Each individually bagged mask and hood has a sterility indicator.

Technical Specification :

Sterile Integrated Hood and Mask



Laboratory Instrument India