



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
50 Coveralls with Hood and Boot

**Product Code :**  
LBNY-00017-11000184



**Description :**

Made with a microporous film that is laminated to a spund bonded polypropylene substrate  
Made from microporous film, garments help protect workers who face on-the-job exposure to a range of non-hazardous aerosols, light liquid splash and dry particles.  
Front zipper closure  
Elastic wrist  
Standard hood with elastic around the face  
Attached skid-resistant boots  
White  
Serged seams  
Thickness: 10 mils.

**Technical Specification :**

50 Coveralls with Hood and Boot



**Laboratory Instrument India**

