



Email : sales@laboratoryinstrumentindia.com

Product Name :
CPE Shoe Covers

Product Code :
LBNY-00017-11000188



Description :

Lightweight all film construction
Fluid resistant
Ideal for short wear applications or transition through cleanroom dress area. Available in blue or white.
Seamless sole with heat-sealed seams
Critical Cover shoe and boot covers slide easily over this shoe cover to speed up donning of final cleanroom cover.
Latex-free.

Technical Specification :

CPE Shoe Covers



Laboratory Instrument India