



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sleeves

Product Code :
LBNY-00017-11000190



Description :

Material is composed of a breathable polyethylene film/spun bond polypropylene nonwoven laminated fabric. Protective sleeves are constructed of a high quality microporous composite material with a low particle count to meet the requirements of controlled environments.

- Static dissipative
- Sonic welded seams
- Heat dissipating material
- Low particle count
- Tapered sleeve with elastic at both ends.

Technical Specification :

Sleeves



Laboratory Instrument India