



Email : sales@laboratoryinstrumentindia.com

Product Name :
Lab Coats

Product Code :
LBNY-00017-11000191



Description :

Strong, lightweight and breathable material
Features serged seams, five-snap closure, elastic wrists, tapered collar, and three pockets
Protection from fluids, light chemical splash, particles and biologicals
Strong, lightweight and breathable material provides protection from fluids, light chemical splash, particles and biologicals without compromising comfort.
Low particle count and high filtration efficiency of 99.9%
Composed of breathable polyethylene film, spunbond polypropylene nonwoven laminated.

Technical Specification :

Lab Coats



Laboratory Instrument India