



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Shoe Covers Serged Seam

**Product Code :**  
LBNY-00017-11000195



**Description :**

Non-conductive polypropylene shoe covers are for use in non-hazardous environments. Ideal for general purpose or for facility visitors.

Serged at the seam

Moderate barrier to fluids, particulates and contamination

Blue.

Available in anti-skid.

**Technical Specification :**

Shoe Covers Serged Seam



**Laboratory Instrument India**