



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Coverall with Elastic Wrist and Ankle and Serged Seams

**Product Code :**  
LBNY-00017-11000204



**Description :**

Protects from grime, dirt and particulates  
Tapered collar neckline and a front zipper for easy donning and doffing  
With elastic wrist and ankle, serged seams.  
Strong, lightweight innovative microporous composite material protects the wearer from grime, dirt, and particulates.  
Pocketless to help control items entering and leaving work area.

**Technical Specification :**

Coverall with Elastic Wrist and Ankle and Serged Seams



**Laboratory Instrument India**