



Email : sales@laboratoryinstrumentindia.com

Product Name :
Critical Cover Coveralls

Product Code :
LBNY-00017-11000206



Description :

Lightweight, durable material lets the wearer stay clean and comfortable all day while providing a limited barrier to particulates and contamination.

Soft polypropylene coveralls are ideal for non-hazardous environments.

Lightweight, durable material

Serged seams

No pockets.

Elastic back

Elastic wrist, ankle and back.

Technical Specification :

Critical Cover Coveralls



Laboratory Instrument India