



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Ankle High Boot Covers

**Product Code :**  
LBNY-00017-11000207



**Description :**

Low particle counts

Boot covers provide anti-skid properties in wet and dry conditions. Two styles to choose from include seamless sonic welded or stress-resistant serged seams.

Available as seamless sole with sonic welded seams or elastic all around with serged seams.

**Technical Specification :**

Ankle High Boot Covers



**Laboratory Instrument India**