



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mens Button-Front Lab Coat

Product Code :
LBNY-00017-11000213



Description :

For ultimate protection from leaks and spills in the lab, look no further than our Lab Coat. Its made out of Performance Blend poplin and comes with three large pockets so that you can keep everything you need close by.

Wash Care: Industrial Laundry - Heavy Soil, Industrial Laundry - Light Soil, Home Wash

Fabric: 5 oz. Performance Blend Cotton

Blend: 80% Polyester/20% Cotton

Primary Closure: Buttons

Pocket: Left interior chest pocket, two lower exterior pockets

No buttons mean fewer repairs and additional protection

Requires minimum ironing thanks to a wrinkle-resistant finish

Interior chest pocket and two lower pockets provide storage for tools and personal items

Finish: Wrinkle Resistant

Closure: Four buttons

Collar: One-piece, notched lapel

Length: 41.5 in. (based on size medium - back)

Gender: Male

Materials: Cotton/Polyester

Side vent openings for easy access to underneath attire

41.5 in. back length on size medium

Fit Details: Fits true to size. Meant to be worn over clothes.

Technical Specification :

Mens Button-Front Lab Coat



Laboratory Instrument India